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Reproducible pain reduction results¹⁻⁹

WOMAC pain subscale change from baseline across studies¹⁻⁹

Pain relief is maintained up to 5 years⁹

In all our studies¹⁻⁹, Arthrosamid® exceeds the MCID** of 9 points¹, and improvement is maintained up to 5 years after single injection⁹.

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
[†] WOMAC or The Western Ontario and McMaster Universities Osteoarthritis Index is a measure of symptoms and physical disability. LS Means are modelled/estimated means. The estimated means are using data from the other visits and also the covariates. ^{**}The minimum clinically important difference (MCID) represents the smallest improvement considered worthwhile by a patient.

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
Referenties: 1. Bliddal, H., et al. (2023) A Prospective Study of Polyacrylamide Hydrogel

Injection for Knee Osteoarthritis: Results From 3 Years After Treatment. *Osteoarthritis and Cartilage*. Vol 31(5): 682-683;

2. Bliddal, H., et al. (2024)

3 year follow-up from a randomized controlled trial of intra-articular polyacrylamide hydrogel injection in subjects with knee osteoarthritis. Poster LB-31, OARSI 2024 World Congress on Osteoarthritis;

3. Neogi, T. (2013) The epidemiology and impact of pain in osteoarthritis. *Osteoarthritis and Cartilage*. Vol 21 (9):

1145-53; 4. Thomson, A., et al. (2021). Synovial

Macrophages in Osteoarthritis: The Key

to Understanding Pathogenesis? *Front. Immunol.* 61

12 :678757; 5. Mathiessen, A., et al. (2017). Synovitis in

osteoarthritis: current understanding with therapeutic

implications. *Arthritis research & therapy*. Vol 19(11):18; 6.

Baker, K., et al. (2010). Relation of synovitis to knee pain

using contrast-enhanced MRIs. *Annals of rheumatic diseases*.

Vol 69(10):1779-83; 7. Data on file; 8. Torres, L., et al.

(2006). The relationship between specific tissue lesions and

pain severity in persons with knee osteoarthritis. *Osteoarthritis and Cartilage*. Vol 14(10): 1033-40; 9.

Christensen, L., et al. (2016). Histological appearance of the

Synovial Membrane after Treatment of Knee Osteoarthritis

with Polyacrylamide Gel Injections: A Case Report. *Journal of*

Arthritis. Vol 5: 217; 10. Bliddal, H., et al. (2021)

Polyacrylamide Hydrogel Injection for Knee Osteoarthritis: A 6

Months Prospective Study. *J Orthop Res Ther*. Vol

6 (2). 1188. ISSN 2575-8241; 11. Arthrosamid®, Instructions

For Use. Release Date March 2022. 10082-003;

12. Bliddal, H., et al. (2022). A Prospective Study of

Polyacrylamide Hydrogel Injection for Knee Osteoarthritis:

Results From 2 Years After Treatment. Poster presented at

OARSI 2022. *Osteoarthritis and Cartilage*. Vol 30(1): S371-

S372; 13. Bliddal, H., et al. (2021). Polyacrylamide

Hydrogel Injection for Knee Osteoarthritis: Results of a 52

Week Prospective Study. *Osteoarthritis*

and Cartilage. Vol 29: S278; 14. Bliddal, H., et al. (2022).

One-year performance of polyacrylamide hydrogel vs.

hyaluronic acid: A randomised controlled study. *Osteoarthritis*

and Cartilage. Vol 30, S1, S371; 15. Christensen, L., et al.

(2016). Synovial incorporation of polyacrylamide hydrogel

after injection into normal and osteoarthritic animal joints.

Osteoarthritis and Cartilage. Vol 24(11):1999-2002; 16.

Grässel, S., et al. (2020). Recent advances

in the treatment of osteoarthritis. *F1000Res*. Vol

9:F1000 Faculty Rev-325; 17. Henriksen, M., et

al. (2018). Intra-articular 2.5% polyacrylamide hydrogel for

the treatment of knee osteoarthritis: an observational proof-of-

concept cohort study. *Clin Exp Rheumatol*. Vol

36(6):1082-85. Epub 2018 Jul 18. PMID: 30148430; 18.

Trigkilidas, D., et al. (2013) The effectiveness of hyaluronic

acid intra-articular injections in managing osteoarthritic knee

pain *Ann R Coll Surg Engl*. Vol 95: 545–551; 19.

Overgaard, A., et al. (2019). Safety of intra-articular

polyacrylamide hydrogel for the treatment of knee

osteoarthritis symptoms: A retrospective case series. *Clin*

Ortho Adv Res. Osteoarthritis and Cartilage. Vol 30(1):

S370-71. 20. Bliddal, H., et al. (2024) A Prospective

Study of Polyacrylamide Hydrogel Injection for Knee

Osteoarthritis: Results From 4 Years After Treatment.

Presented at EORS 2024, Aalborg and Orthopaedic

Proceedings of the Bone & Joint Journal. 21. 10-Year

Safety Follow-Up of Intra-Articular 2.5%

Polyacrylamide Hydrogel (iPAAG) in Knee

Osteoarthritis (OA) Presented at: WCO-IOF-ESCEO

2025, Rome; H. Bliddal, et al. 21. A Mechanistic and

Clinical study of intra-articular Arthrosamid for knee

osteoarthritis Presented at the ISAKOS congress

2025 ; Sharon Owen, et al. 22. Polyacrylamide

Hydrogel for Knee Osteoarthritis: 5-Year Results from

a Prospective Study Presented at WCO-IOF-ESCEO

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ARTHROSAMID[®]
by
contura
ORTHOPAEDICS